# DentaPic SONIC

## SCHALL-ZAHNREINIGER

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**EN**

## **IMPORTANT SAFEGUARDS**

Please read the safety instructions before using this product.  
The safety instructions provided here are designed to help

you use this product safely and correctly, and also to prevent you from injuring yourself or others. When using electrical products especially when children are present, basic safety precautions should always be followed, including the following.

## **WARNING**

- PLEASE READ ALL INSTRUCTIONS BEFORE USING THIS APPLIANCE AND RETAIN THESE INSTRUCTIONS FOR REFERENCE.
- This device may be used by children from 8 years of age, or by persons with reduced physical, sensory or mental capabilities or lack of experience and/or knowledge, only if they have been supervised or instructed in the safe use of the device and have understood the hazards arising therefrom. Children should not play with the device. Cleaning and user maintenance must not be performed by children unless they are over 8 years old and supervised. The device is to be kept away from children younger than 8 years of age.
- This product is solely intended for cleaning teeth. Do not use any attachments other than those recommended by the manufacturer.
- If you experience discomfort or pain in your mouth, cease use of the device and consult a dentist.
- **IMPORTANT:** The device is not a substitute for daily

brushing.

- If you have had oral or gum surgery in the past 2 months, please consult your dentist before using this product.
- If you experience heavy bleeding after using this device, or if you still bleed after using it for a week, cease use immediately and consult your dentist.
- Consult your dentist if you wear a pacemaker or other medical implant.
- Do not clean the product in a dishwasher.
- This product is for indoor use only, and is not intended for outdoor use.
- If not used for a longer period of time, remove the batteries to avoid danger or damage to the device, such as that caused by leakage of batteries.
- Make sure that the skin, eyes and mucous membranes never come into contact with the battery fluid/acid. In case of contact with the battery fluid, immediately rinse the affected areas with plenty of clean water and consult a doctor.
- Consult your doctor before using the DentaPic Sonic if you have any medical issues.

## **INSTRUCTION MANUAL**

- Can gently remove stains, plaque and tartar.
- Can clean teeth and interdental spaces.
- Can promote the health of gums.

**To install the battery:**

- Locate the battery compartment cover on the lower side of the device.
- Turn and pull the battery compartment cover off.
- Insert 1x AAA battery (not included) and close the battery cover again by pressing and turning until the two arrows are aligned.

**Application:**

- Select the attachment you wish to use; either the stainless steel attachment or the silicone attachment.
- To change the accessory, simply pull up the attachment in use, remove it, and replace with the desired attachment.
- Turn on the DentaPic Sonic by pressing the button located at the side twice. Press once to turn the lights on only. Then clean your teeth and interdental spaces by simply gently pressing down on to stains, plaque and interdental spaces and holding in place!
- Rinse with water after use!

**Please note the following important tips for best results:**

- To remove tartar, select the stainless steel attachment. Application should not be performed more than once per week.
- To remove plaque and discolourations, select the silicone attachment. Application can be as needed (even several times per day).
- Never apply the device directly to the tooth, but instead at a slight angle. Hold the tip on the area to be cleaned as shown in Illustration 1.

**IMPORTANT!** To avoid injury during use, perform all movements slowly and do not apply pressure. Gently apply the tip to the area to be cleaned. Clean only along the gum line - do not apply directly to the gums.

**Note for crown, veneer and brace wearers:**

Consult your dentist prior to use on (especially older) crowns, veneers and braces.

**Caution: The device must NOT be used with dental implants.**

**Cleaning:**

**IMPORTANT:** *Do not immerse in water.*

- Make sure that the DentaPic Sonic is turned off before you begin to clean it.

- Clean the accessories by removing them from the device and washing with water or cleaning with alcohol swabs.

**Technical data:**

Battery: 1x 1,5 V AAA; Power: 1,5 V  2 W

(battery not included)

**IMPORTANT INFO:** DentaPic Sonic does NOT replace brushing, but instead serves as an additional dental care product. Consult your doctor before using the DentaPic Sonic if you have any medical concerns.

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**EN:** Do not dispose this product at the end of its service life through normal domestic waste. It has to be disposed through a collection agency that recycles electric and electronic devices. This is pointed out by the symbol on the product, in the instruction manual, and on the packaging. Please inquire about local collection agencies that may be operated by your distributor or your municipal administration. By recycling, making use of the materials, or other forms of recycling old units you are making an important contribution to protecting our environment.

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**CH:** MediaShop AG | Industriering 3 | 9491 Ruggell | Liechtenstein
**EU:** MediaShop GmbH | Schwarzottstraße 2a | 2620 Neunkirchen | Austria